

Target Sequence Cloning Protocol

(Standard de-salted oligos are sufficient)

PX330-based plasmids, including PX458-462 – SpCas9 (or SpCas9n D10A nickase) + single guide RNA:

To clone the guide sequence into the sgRNA scaffold, synthesize two oligos of the form:

```

5' - CACCGNNNNNNNNNNNNNNNNNNNN - 3'
3' -      CNNNNNNNNNNNNNNNNNNNCAA - 5'

```

PX260 and PX334 – SpCas9 (or SpCas9n D10A nickase) + CRISPR array + tracrRNA:

To clone the guide sequence into the sgRNA scaffold, synthesize two oligos of the form:

```

5' - AACNNNNNNNNNNNNNNNNNNNNNNNNNGT - 3'
3' -      NNNNNNNNNNNNNNNNNNNNNNNNCAAAT - 5'

```

* * * * *

Oligo annealing and cloning into backbone vectors:

1. Digest 1ug of plasmid with *Bbsl* for 30 min at 37°C:

1 ug	Plasmid
1 ul	FastDigest <i>Bbsl</i> (Fermentas)
1 ul	FastAP (Fermentas)
2 ul	10X FastDigest Buffer
X ul	ddH ₂ O
<hr/>	
20 ul	total

2. Gel purify digested plasmid using QIAquick Gel Extraction Kit and elute in EB.

3. Phosphorylate and anneal each pair of oligos:

1 ul	oligo 1 (100uM)
1 ul	oligo 2 (100uM)
1 ul	10X T4 Ligation Buffer (NEB)
6.5 ul	ddH ₂ O
0.5 ul	T4 PNK (NEB)
<hr/>	
10 ul	total

Anneal in a thermocycler using the following parameters:

37°C	30 min
95°C	5 min and then ramp down to 25°C at 5°C/min

4. Set up ligation reaction and incubate at room temperature for 10 min:

X ul	<i>Bbsl</i> digested plasmid from step 2 (50ng)
1 ul	phosphorylated and annealed oligo duplex from step 3 (1:200 dilution)
5 ul	2X Quickligation Buffer (NEB)
X ul	ddH ₂ O
<hr/>	
10 ul	subtotal
1 ul	Quick Ligase (NEB)
<hr/>	
11 ul	total

5. (optional) Treat ligation reaction with PlasmidSafe exonuclease to prevent unwanted recombination products:

11 ul	ligation reaction from step 4
1.5 ul	10X PlasmidSafe Buffer
1.5 ul	10mM ATP
1 ul	ddH ₂ O
<hr/>	
15 ul	total

Incubate reaction at 37C for 30 min.

6. Transformation